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Temperature control jacket

ZTJ 3020

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Consistent temperature during viscosity determination ensures meaningfully results.

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Application areas

- for paint and varnish manufacturers, its customers and test laboratories
- for paints, varnishes and similar liquids
- the test can be carried out in the laboratory and during production
- applicable with all common flow cups corresponding to the standards

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Standard delivery

- 1 temperature control jacket for flow cups
- 1 closed double featured vessel for indirect tempering
- 2 tube accesses to temper liquids
- 1 spirit level for levelling
- 1 certificate of manufacturer

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Options

- ZFC 3010-3011 as well as ZFC 3013 Flow cups
- ACC050 insertion thermometer up to 400°C (752°F)

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Handling

- thanks to the adjustable feet and spirit level, the temperature control jacket is easily to levelling
- put the chosen flow cup into the temperature control jacket
- for tempering you need additionally a circulation thermostat
- connect the circulation thermostat by use of two heating/cooling tubes on the tube access of the temperature control jacket so that you have a closed circuit
- adjust the required temperature at the circulation thermostat and begin with tempering
- start the viscosity determination at the required temperature in accordance to the separate instruction of the flow cup and the respecting standard
- after finishing viscosity determination, firstly remove one heating/cooling tube, cover the tube access with your thumb, remove the second heating/cooling tube and pour the temperature control liquid in the circulation thermostat

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Technical specification

version	application / standards	flow cup
3020.ISZ	for cylindrical Zehntner-Flow cups according to DIN 53211 or ASTM D 1200 with adaptor for flow cups according to ISO 2431	
3020.DFK	for conical Zehntner-Flow cups according to DIN 53211 or ASTM D 1200	
3020.FDF	adaptors for flow cups of other manufactures according to DIN 53211 or ASTM D 1200	

material: anodised aluminium
 dimensions (LxWxH): 180 mm x 210 mm x 260 mm (7.1" x 8.3" x 10.2")
 weight: 3020.DFK: 750 g (1.65 lbs), 3020.ISZ: 800 g (1.76 lbs), 3020.FDF: depending on adaptor; approx. 800 g (1.76 lbs)
 warranty: 2 years